



Ready Matic

Digital Blood Pressure Monitor

CH-484

Jumbo Digits Display for Easy Reading
Ultra-large memory capacity 100 memories
offering data for you or reference of your doctor.

Unique heart display
The unique lively figure
design is based on WHO
Blood Pressure Classification
Chart for fast and easy reading.



Twin Memory Zones
Total 100 memories

- Fuzzy Logic
- Measuring method: Oscillometric
- With big LCD display
- Jumbo digits display for easy reading
- One touch: easy to be operated
- Big memory space: 2 x 50 readings
- Intelligent error detection and display
- With low battery display
- Unique heart display to reflect the measuring condition immediately
- Automatic shutdown device

| Figure | Range | Systolic | Diastolic |
|--------|--------------|----------------|--------------|
| | Hypertension | > 140 mmHg | > 90 mmHg |
| | High-Normal | 130 ~ 140 mmHg | 85 ~ 90 mmHg |
| | Normal | < 129 mmHg | < 84 mmHg |
| | Hypertension | < 90 mmHg | < 50 mmHg |

| | |
|-------------------------|--|
| Measuring method | Oscillometric method |
| Measuring range | Pressure: 20-280 mmHg Pulse rate: 40-200 pulses/min |
| Accuracy | Pressure: ±3 mmHg Pulse: ±5% of reading value |
| Memory capacity | 2x50 readings |
| Pressure inflation | Automatic pump actuated device |
| Pressure deflation | Automatic linear pressure deflation valve |
| Pressure release device | Automatic solenoid venting valve |
| Sensor | Semiconductor pressure sensor |

Specifications

| | |
|----------------------------------|--|
| Electric power | DC6V (four LR06 "AA" alkaline batteries) |
| Operation temperature / humidity | 10 - 40°C (50 - 140 F) 30% - 85% R.H. |
| Storage temperature / humidity | - 5 - 50°C (23 - 122 F) 15% - 90% R.H. |
| Weight | Approx. 442 grams (without batteries) |
| Arm cuff. | Circumference 23~33 cm (9"~13") |
| Dimensions | 135 (L) x 115 (W) x 75 (H) mm |
| Accessories | Storage bag, user manual, AA Alkaline battery x 4 pcs; Cuff. |

NOTE : Specifications are subject to change without further notice.

Your Area Distributor



Exclusively Imported By :

Mehar Healthcare Corporation

Admn. Off :16/89, Lane No. 4, Faiz Road, Karol Bagh, New Delhi - 110005
Phone : 011-23520753, 23539529, 23864115, 23861140
Fax : 011-25812657 Mobile : 98101-00623
E-mail : meharhealthcare@rediffmail.com
www.meharhealthcare.com